

G3

Perfect for perfectionism

Quality: Over 200KG weight supplied gear motor with big torque force , long lifespan and steady quality

Profession: Top-mounted tray with intelligent controls

Affordable Luxury: Six led bulbs lamp dual-color

Humanization: Humanized functions for high effective work experience

Fashion: Fashionable designed side box which meets the high end taste



*The G3 name is correspond to HJ-668A in the registration certificate

The elegant side box adds flexibility when choosing a treatment maneuver during a day



Upper hanging instrument tray, light and flexible pull bar make the operation more convenient



Cup holder and tissue box




Side box circuit optimization, separation of water and electricity, easy to maintain



The Armrest can be rotated and convenient for patient

Seat Size	Large Width - 54cm , Seat - 64cm Backrest
Package	143*102*115cm , G.W : 255KG
Container Volume	16units/20GP , 32units/40HQ

Major Parts
1. Water heating system
2. Spittoon flushing and cup filling
3. Quiet electric lift chair
4. Fashionable hanging main table
5. Chair movement system with 3 memory positions
6. LED film viewer
7. Assistant control panel
8. Rotatable ceramic spittoon
9. Shadowless LED sensor dental light
10. Multifunction pedal switch
11. Pure water bottles, 600ml
12. Glass tray and tissue box
13. Dentist stool with backrest
14. Imported polyether inner tubing



Technical Specification
Supply voltage : AC110V/220V 60Hz/50Hz
Motor voltage : DC24V
Air pressure : 0.5MPa - 0.8MPa
Water pressure : 0.2MPa - 0.4MPa

Dentist Stool



01



02



03

Optional Dental Lamp



L1



L2



L3

Optional Accessory



Handpiece Set



Built-in Scaler



Built-in LED
Curing Light



Intra-oral
Camera Set



Medical
Compressor



Top-mounted
Main Table

Upholstery Color

Microfiber Upholstery



PU Upholstery

