

G7

Imported motor
9 memory system
Extra stainless steel tray
Injected Plastic
6 LED Dual Color sensor lamp
Laser printed touch button assistant tray



**EXQUISITE LUXURY
DESIGN MAKES
THE CHAIR AS PREMIUM
PRODUCT**



Side box circuit optimization,
separation of water and electricity,
easy to maintain



The Spittoon can be rotated and
remove for easy cleaning

Seat Size	Large Width - 54cm , Seat - 64cm Backrest
Package	143*102*115cm , G.W : 255KG
Container Volume	16units/20GP , 32units/40HQ

Major Parts

1. Water heating system
2. Spittoon flushing and cup filling
3. Quiet electric lift chair(supported by Denmark motors)
4. Fashionable control panel
5. Chair movement system with 9 memory positions
6. LED film viewer
7. Assistant control panel
8. Rotatable ceramic spittoon
9. Shadowless LED sensor dental light
10. Multifunction pedal switch
11. Pure water bottles 600ml
12. Dentist stool with backrest
13. Imported polyether inner tubing



Technical Specification

Supply voltage : AC110V/220V 60Hz/50Hz

Motor voltage : DC24V

Air pressure : 0.5MPa - 0.8MPa

Water pressure : 0.2MPa - 0.4MPa

Dentist Stool



01



02



03

Optional Dental Lamp



L1



L2



L3

Optional Accessory



Handpiece Set



Built-in Scaler



Built-in LED
Curing Light



Intra-oral
Camera Set



Medical
Compressor



Top-mounted
Main Table

Upholstery Color

Microfiber Upholstery



PU Upholstery

