

H5

Perfect for perfectionism

Quality: Over 200KG weight supplied

gear motor with big torque force , long

lifespan and steady quality

Humanization: Humanized functions for

high effective work experience

Fashion: Fashionable designed side box

which meets the high end taste

Imported inner tubing



Adjustable shadowless dual color dental lamp



The 90 degree rotation spittoon is convenient for the patient to gargle, smooth surface and easy to clean



Due to the incorporation of ergonomic principles into the design, the work has become very simple and intuitive, and the burden on the body has been minimized



Integrated suction filter, easy for cleaning



A TRUSTWORTHY CLASSIC

Seat Size	Normal Width - 50cm , Seat - 57cm Backrest
Package	135*93*108cm , G.W : 190KG
Container Volume	20units/20GP , 40units/40HQ

Major Parts

1. Water heating system
2. Spittoon flushing and cup filling
3. Quiet electric lift chair
4. Fashionable main table
5. LED film viewer
6. Rotatable ceramic spittoon
7. Assistant control panel
8. Shadowless LED sensor dental light
9. Multifunction pedal switch
10. Pure water bottles, 600ml
11. Dentist stool with backrest
12. Imported polyether inner tubing



Technical Specification

Supply voltage	: AC110V/220V 60Hz/50Hz
Motor voltage	: DC24V
Air pressure	: 0.5MPa - 0.8MPa
Water pressure	: 0.2MPa - 0.4MPa

Dentist Stool



01



02



03

Optional Dental Lamp



L1



L2



L3

Optional Accessory



Handpiece Set



Built-in Scaler



Built-in LED
Curing Light



Intra-oral
Camera Set



Medical
Compressor



Top-mounted
Main Table

Upholstery Color

Microfiber Upholstery



PU Upholstery

