

P2

Perfect Solution For High Demand

Firmness: Avant-garde designed side box

Comfort: Grand & ergonomic designed seat

Excellence: Memory system

Profession: LED Sensor dual-colored lamp

Quality: Over 200KG weight supplied

gear motor with big torque force , long

lifespan and steady quality

Imported inner tubing



Side box with 90 degrees rotation, the best solution for 2-hand or 4-hand treatment



Rotatable touch button panel

SUPREME TREATMENT CREATING THE PERFECT EXPERIENCE



Seat Size	Large Width - 54cm , Seat - 64cm Backrest
Package	143*102*115cm , G.W : 290KG
Container Volume	10units/20GP , 32units/40HQ

Major Parts
1. Water heating system
2. Spittoon flushing and cup filling (supported by Italian solenoid valves)
3. Quiet electric lift chair (supported by Denmark motors)
4. LCD display control panel
5. Chair movement system with 9 memory positions
6. Metal backrest seat frame
7. Assistant control panel
8. Rotatable ceramic spittoon
9. Shadowless LED sensor dental light
10. Multifunction pedal switch
11. Pure water bottles 600ml
12. Glass tray and built-in tissue box
13. Deluxe dentist stool
14. Imported polyether inner tubing
15. Side box unit 90°rotatable



Technical Specification
Supply voltage : AC110V/220V 60Hz/50Hz
Motor voltage : DC24V
Air pressure : 0.5MPa - 0.8MPa
Water pressure : 0.2MPa - 0.4MPa

Dentist Stool

01

02

03

Optional Dental Lamp

L1

L2

L3

Optional Accessory

Handpiece Set

Built-in Scaler

Built-in LED Curing Light

Intra-oral Camera Set

Medical Compressor

Top-mounted Main Table

Upholstery Color

Microfiber Upholstery

PU Upholstery