

QL2028 IV

牙科综合治疗机 / DENTAL UNIT

完美于形

Perfect in form





功能介绍 / FEATURES



转动机箱
Rotary side box


旋转灵活，增大助手位的活动空间，便于副手操作。

Easy to move, more space for the assistant, make the assistant operation easier.



副控台
Assistant control

八键多功能操作面板；
三用枪（热）；
可调节强吸接头；
Multi- function control panel with 8 buttons;
3-way syringe (heat);
Adjustable strong suction;



外置可拆卸式吸唾过滤器
Built out and detachable
suction filter

快速拆装，易于清洁，安装高质量强弱吸管。

Easy to clean with good quality strong and weak suction tube.



全新上挂器械盘
Top-mounted delivery tray

精心设计的上挂式器械盘，使用方便，操作轻松自如，使用过程几乎没有后拽感。

Well-designed top-mounted tray, easy to operate.



功能性的控制面板
Full function control panel

LCD 液晶显示屏显示工作状态，3 位医生 9 组记忆位设置键，吐痰联动椅位快捷键，口腔灯快捷键，观片灯快捷键，术后复位快捷键。

With LCD screen displaying the working status, 3 doctors with 9 sets of chair position memory, spitting position setting, operation light control, x-ray film viewer control, zero position control.



舒适联动治疗椅设计及机椅互锁系统
Comfortable patient chair with chair
interlock system

人体工程学设计，联动椅架，高档超纤皮垫，配机椅互锁系统，更舒适，更安全。

Synchronized chair with big and handmade sewed ultra-fiber leather upholstery, with chair interlock system, more safe and comfortable.



高档旋转式玻璃痰盂
Rotary glass cuspidor

玻璃痰盂质感强，易拆卸，清洗方便，吐痰舒适，独特的水流疏导设计，有效防止水流飞溅。

Rotary and detachable spittoon,easy to clean.



感应 LED 柔和双色灯
LED sensor light with yellow and
white light mix

黄光、白光、柔和混合光三种光源模式自由切换，光斑柔和清晰均匀，感应手动两用开关，亮度强至弱可调节，手柄可拆卸，光源分为两组黄光以及四组白光。

More tender light: 6 light source with 4 white and 2 yellow; Easily switch to yellow light mode, white light mode and mix light mode; Super clear light spot, brightness adjustable from strong to weak; Detachable and sterilizable handle.



豪华医生椅
Luxury dentist stool

宽大坐垫以及椅背，更加舒适，更加安全。

Big cushion with strong metal frame, more safe and comfortable.