

QL2028I (2024)

Stylish appearance
Stable performance



Easy to control
Clear functions at a glance

Strictly control every detail

Scientifically designed backrest support

New PU foam armrests with delicate touch, making it easy to experience high-end treatment. Curved design provides comfortable support to the neck and prevents the head from shaking left and right. Comfortable backrest, not tired after lying down for a long time, allowing patients to experience comfortable treatment, carefully designed to protect the patient's waist. Lowest chair position is 40cm, making it convenient for the elderly and children to get in and out of the chair.

Rotating glass cuspidor

Can be rotated 180° to get closer to the patient. Deeper design of the cuspidor effectively prevents water splashing. Spittoon is removable, making cleaning and disinfection more convenient.

External suction

Twist-type compression sealing structure, making it easier for operation. Filter is on the outside of the side box and can be seen visually, making it easier to clean the daily filter.

Independent pure water supply system

Provide pure water supply to the instrument tray 3-way syringe and high and low speed handpieces.



8-key assistant control

Cuspidor flushing, water supply for cup, water heating, spitting position shortcut keys.

Strong and weak suction system, 3-ways syringe (warm).

High quality sensor dental light

6 high color rendering index imported LED lamp beads.
3 light source modes to meet various needs.
2 operating modes, sensor/touch.
8-25k Lux brightness adjustable.
100,000 times life test.

light spot is uniform and non-dazzling; constant current design, no flicker and more eye-friendly.
Patented technology, adjustable elasticity of movable joints.



Colors of splendor
Choose at will

PU cushion

P01

Light Blue

P02

Apple Green

P03

Light Green

P04

Diamond Blue

P05

Yellow

P06

Orange

P07

Red

P08

Dark Purple

P09

Purplish Red

P10

Dark Blue

P16

Seaweed Green

Microfiber Cushion

FM01

Beige

FM02

Red

FM03

Blue

FM04

Coffee

FM05

Light Green

FM06

Dark Blue

FM07

Apple Green

FM08

Orange

Installation diagram

