

QL2028I (2024)

Stylish appearance
Stable performance



Easy to control
Clear functions at a glance

Strictly control every detail

Scientifically designed backrest support

New PU foam armrests with delicate touch, making it easy to experience high-end treatment. Curved design provides comfortable support to the neck and prevents the head from shaking left and right. Comfortable backrest, not tired after lying down for a long time, allowing patients to experience comfortable treatment, carefully designed to protect the patient's waist. Lowest chair position is 40cm, making it convenient for the elderly and children to get in and out of the chair.

Rotating glass cuspidor

Can be rotated 180° to get closer to the patient. Deeper design of the cuspidor effectively prevents water splashing. Spittoon is removable, making cleaning and disinfection more convenient.

External suction

Twist-type compression sealing structure, making it easier for operation. Filter is on the outside of the side box and can be seen visually, making it easier to clean the daily filter.

Independent pure water supply system

Provide pure water supply to the instrument tray 3-way syringe and high and low speed handpieces.



8-key assistant control

Cuspidor flushing, water supply for cup, water heating, spitting position shortcut keys.

Strong and weak suction system, 3-ways syringe (warm).

High quality sensor dental light

6 high color rendering index imported LED lamp beads.
3 light source modes to meet various needs.
2 operating modes, sensor/touch.
8-25k Lux brightness adjustable.
100,000 times life test.

light spot is uniform and non-dazzling; constant current design, no flicker and more eye-friendly.
Patented technology, adjustable elasticity of movable joints.



Colors of splendor
Choose at will

PU cushion



P01
Light Blue



P02
Apple Green



P03
Light Green



P04
Diamond Blue



P05
Yellow



P06
Orange



P07
Red



P08
Dark Purple



P09
Purplish Red



P10
Dark Blue



P16
Seaweed Green

Microfiber Cushion



FM01
Beige



FM02
Red



FM03
Blue



FM04
Coffee



FM05
Light Green



FM06
Dark Blue



FM07
Apple Green



FM08
Orange

Installation diagram

