

S610

DENTAL UNIT


Simple and comfortable
Reliable to use



Keep Improving Casting Quality

- ▶ **Assistant rack silicone pad**
Anti scratch, easy to clean and disinfectable

- ▶ **Instrument tray silicone pad**
Anti scratch, easy to clean and disinfectable



- ▶ **The chassis can rotate 45°**
Convenient four handed operation



- ▶ **Aluminum alloy instrument tray (interlocking)**
The handle can be disassembled for disinfection, and the hanging rack can be rotated



- ▶ **LED oral light**
Anti resin curing mode, induction on/off brightness adjustment



- ▶ **External water storage bottle device**
One click conversion/switch, convenient and fast




- ▶ **Detachable strong weak suction filtration device**
Easy to clean

S610 DENTAL UNIT

Engineering electric chair Comfortable seat cushion design

The seat design is based on the human body shape, providing perfect comfort and support, especially for patients undergoing long-term treatment.

Keep Improving Casting Quality



- **Double joint headrest**
Folding design makes patients more comfortable and at ease during treatment




- **Rotable armrest design**
Aluminum alloy material / engineering plastic cover



- **Emergency stop**
One click stop, safe, and convenient



- **The chasis can rotate 45°**
Convenient four handed operation



- **Silent electric chair(Protect)**
Designed with a minimum seating position of 380mm, convenient for children and elderly patients



- **Ceramic spittoon, automatic water supply system**
Detachable/rotating, flushing and mouthwash quantitative water supply



- **Multifunctional foot pedal**
Control up and down/supine/ cupping/ water supply/ oral light / spitting position/ mobile phone working water status/ speed



- **Water/gas/electricity main switch device**
One click on/off, convenient and fast (anti suction device)

