

S620 DENTAL UNIT

Simple and comfortable
Reliable to use

Striving for excellence, casting quality



LED oral light

Removable handle for cleaning and disinfection



LED oral light(Digital brightness/ color temperature)

Induction on/off brightness adjustment (manual switching of three light source modes: white light, yellow light, and yellow white mixed) to reduce eye fatigue and achieve accurate color comparison



Assistant rack combination device(Delayed shutdown)

Aluminum alloy handle detachable for cleaning and disinfection



Assistant rack silicone pad

Anti scratch, easy to clean, and disinfectable




External water storage bottle device

One click conversion/switch, convenient and fast



Detachable strong weak suction filtration device

Easy to clean



Instrument tray silicone pad

Anti scratch, easy to clean, and disinfectable

Aluminum alloy instrument tray(interlocking)

The handle can be disassembled for disinfection, and the hanging rack can be rotated

Keep Improving Casting Quality

The light pulling tool is adjustable in position and can be positioned and locked, which conforms to ergonomic principles. It is easy to operate and reduces wrist fatigue.

S620

DENTAL UNIT

Engineering electric chair Comfortable seat cushion design

The seat design is based on the human body shape, providing perfect comfort and support, especially for patients undergoing long-term treatment.



Keep Improving
Casting Quality

- **Double joint headrest**
Folding design makes patients more comfortable and at ease during treatment
- **Rotable armrest design**
Aluminum alloy material/ engineering plastic cover

Striving for excellence , casting quality



- **Glass spittoon, automatic water supply system**
Detachable/rotating, flushing and mouthwash quantitative water supply



- **The chassis can rotate 45°**
Convenient four handed operation



- **Multifunctional foot pedal**
Control up and down/ supine/ cupping/ water supply/ oral light/ spitting position/ mobile phone working water status/speed



- **Water/gas/electricity main switch device**
One click on/off, convenient and fast (anti suction device)



- **Silent electric chair (Protect)**
Designed with a minimum seating position of 380mm, convenient for children and elderly
- **Emergency stop**
One click stop, safe, and convenient

