

QL2028II

臻享款



稳 | 从容所托

3 点式椅架结构 | 350kg 极限承重 | 精工低噪电机 | 5mm 加厚一体成型金属靠背

舒 | 无感亲肤

微感启停系统，臻享无感缓启缓停
曲面头枕 + 超宽蝶形人体工学设计
超纤抗菌亲肤皮料，密合精缝工艺

亮 | 8珠进口LED双色口腔灯

超越明亮 | 三重光源 | 感应 / 手控双模式
4000k-5500k 可调色温
8000-45000LUK 明亮照度
白光 / 黄光 / 混光，一键切换

巧 | 灵动设计

更安全：遇阻回弹 + 机椅互锁 + 封闭电源模块
更灵动：45°可旋转机箱 + 半自动消毒 + 快拆设计
更可靠：进口水气管路 & 电磁阀
更轻量：自研多功能脚踏，精准减轻脚部负担

效 | 高效协作

685mm×430mm 全景工作台，搭载 17 键按键式控制面板及 9 组记忆椅位（全新角度传感技术）。9 键多功能联动助手台，6 位可旋转灵活挂架。更有椅架联动吐痰功能。可调节强弱吸手柄，支持 4 手操作。



QL2028II PRIME



Smart. Flexible Design.

Imported Air & Water Tubing and Solenoid Valves

Multi-Function Foot Pedal for easy treatment

Chair—Unit Interlock

Sealed Power Module

Backrest Anti-collision

Semi—Auto Disinfection

45° Rotatable Side Box

Quick—Release Design

Stable. Trusted. Reliable.

3-Point Frame Structure | Up to 350 kg Load Capacity | Quiet Precision Motor | 5 mm One-Piece Metal Backrest

Comfort. Skin-Friendly Feel.

Soft Start & Stop System for Smooth, Gentle Motion

Curved Headrest + Wide Ergonomic Design

Skin-Friendly Antibacterial Upholstery

Efficient. Collaborative System.

Panoramic Work Surface (685 × 430 mm), carrying 17-Key Control Panel and 9 Memory Chair Positions (Advanced Angle Sensor Technology). 9-Key Assistant Control Panel (Specifications subject to actual product). 6-Position Rotatable Holders.

* Spittoon Rinse Key

* Adjustable Suction Handles

* Supports Four-Handed Dentistry

Bright. 8-LED Two-Color Dental Light.

Adjustable Color Temperature: 4000K—5500K

Illuminance: 8,000—45,000 Lux

White / Yellow / Mixed Light — One-Touch Switching

Exceptional Brightness

Triple Light Source

Sensor & Touch Control

