

Customized Clinic Opening Plan:



Dental Unit Features:

- 1. Sensor operating light
- 2. Water heating system
- 3. Auto spittoon flush and cup filer sensor control system
- 4. Assisitant control box with 2 tubes
- 5. 3 way syringe(cold/hot water)2 pcs
- 6. Saliva ejector with high and low suction
- 7. 24Vnoiseless DC motor chair with one arm
- 8. Instrument tray with air brake with 4 tubes
- 9. Doctor chair
- 10. Multifunction foot pedal switch
- 11.X-film viewer
- 12.600ML Water tank(air pressurized)(built out)

The Clinic opening package included:

- 1.Economical dental unit A1.
- 2. 17 inch monitor with intraoral camera.
- 3. 32Leconomical air compressor.
- 4. Built-in ultrasonic scaler.
- 5.Built-in curing light.
- 6.2pcs high speed handpiece and 1 set of low speed handpiece.
- 7.Clamp on dental chair type teeth whitening unit.

TJ2688 A1

牙科综合治疗台
CONTROLLED INTEGRAL DENTAL UNIT





简约，素雅，经济实惠

SIMPLE, ELEGANT, ECONOMICAL

大尺寸器械盘，配有防污、防尘塑料垫，方便术后清洁和护理。
LARGE SIZE INSTRUMENT PLATE, EQUIPPED WITH ANTIFOULING, DUSTPROOF PLASTIC PAD, EASY TO CLEAN AND CARE AFTER OPERATION.



全方位LED口腔灯
OMNI-DIRECTIONAL LED ORAL LAMP



易清洗痰盂
EASILY CLEAN SPITTOON


一体式箱体和痰盂设计，更加简约实用。
Integrated box and spittoon design, more simple and practical.

选材优良，模具化生产，性能可靠，操作方便，是经济实惠的牙科设备。
Excellent material selection, mold production, reliable performance, convenient operation and economical dental equipment.

标准配置

THE STANDARD CONFIGURATION

出厂配置 STANDARD CONFIGURATION	TJ2688 A1
提拉式扶手 LIFT-DRAW STYLE ARMREST	
可旋转豪华扶手 ROTATABLE LUXURY ARMREST	
可拆卸舒适款扶手 DETACHABLE COMFORTABLE ARMREST	☐
全电脑控制，静音低压直流电机驱动 COMPUTER CONTROLLED NOISELESS LOW VOLTAGE DIRECT CURRENT MOTOR	☐
冲痰，漱口定量给水自动控制系统 MOUTHWASH WATER SUPPLY FLUSHING SPITTOON AUTOMATIC CONTROL SYSTEM	☐
椅位记忆功能 MEMORY FUNCTIONS	
三用喷枪（一热一冷）2支 3 WAY SYRINGE (COLD / HOT)	☐
全方位LED口腔灯强弱两档可感应，可手动两用 STRONG AND WEAK TWO SET LED SENSOR OPERATION LAMP	☐
LED观片灯 LED DENTAL FILM VIEWER	☐
可拆洗玻璃痰盂 DETACHABLE GLASS CUSPIDOR	
辅助控制系统 ASSISTANT CONTROL SYSTEM	
强，弱吸唾装置 STRONG AND WEAK SUCTION	☐
多功能脚踏 MULTI-FUNCTIONAL FOOT PEDAL	☐
圆形脚踏 CIRCULAR FOOT PEDAL	
医生椅 DENTIST STOOL	☐
进口水气管 IMPORTED WATER PIPE	
豪华超纤皮 LUXURY MICROFIBER LEATHER	
旋转侧箱 TURNTABLE SIDE BOX	
可左手操作 LIFT-DRAW STYLE ARMREST	



TJHS-1



Specification:

Bearable weight:350kg
Stand lifting range:12cm Stand
height:42-54cm
Product size:56*56*66cm Packing size:59*
58*26cm Gross weight:10.5kg
Upholstery:PU Leather

Features:

1. Special Designed for dentists.
2. Relieve lumbar muscle strain.
3. Reduce fatigue and nerve compression.
4. Prevention of frozen shoulder and carpal tunnel or vascular syndrome.
5. Solve the thigh blood flow is poor.

TJ-EQ6 Intraoral Camera

17 inch HD intraoral
camera with WiFi

- 1. Smart screen projection with mobile phone
- 2. USB, HDMI output
- 3. Wifi connects to the computer

Packing size: 54*44*17cm Gross weight: 6kg





Built in scaler

Features:

1. Special design for dental unit, built in type.
2. Function: Scaling, Perio.
3. Detachable handpiece with light, better visibility, and autoclavable.
4. Handpiece illuminance 8000 lux
5. Automatic frequency tracking, digitally controlled
6. Shell of the handpiece is made of special imported plastic from America, not easy to be cracking or deformation under pressure.
- 7.5 tips attached: G1x2, G2, G4, P1

Curing light

Built-in Connect to Dental Chair CuringLight

Features:

1. Curing with broad-spectrum blue-violet light (385-515nm), it is possible to cure all composite materials and adhesives.
2. Connect to dental chair, with tube, no need charging.



KY-1EW

Air compression

Feature:

Voltage:	220v50hz
Electric current:	4.8APower:545w
Frequency:	70Lmin
Start pressure:	0.5Mpa
Capacity pressure:	0.8Mpa
Air tank volume:	32L
Noise:	55-62dB
G.W.:	27KG
Packing size:	43*43*70cm

ONE TO ONE AIR COMPRESSOR



Air.com

STU-5 HANDPIECE SET

(2high speed handpieces & llow speed handpiece)



High speed handpiece

Air pressure :	0.2~0.22Mpa(2holes) 0.25~0.27Mpa(4holes)
Rotation speed :	More then 300.000rpm
Chuck type :	Push button
Bur application :	1.9mm~1.6mmx21mm~ 23mm(diameterxlength)
Noise :	Less then 68dB
Spray :	1 way spray



Low speed contra angel

Air oressure :	0.3Mpa
Rotation speed :	About 20,000rpm
Noise :	≤70d
Bchuck type :	Push type
Bur applicaation :	2.35mm
Spray :	External channel



Low Speed Handpiece

Air oressure :	0.3Mpa
Rotation speed :	About 20,000rpm
Noise :	50dB
Ratio :	1:1
Chuck type :	Wrench
Bur application :	2.35mm
Spray :	External channel



Low Speed Air Micro Motor

Air pressure :	0.3Mpa
Rotation speed :	About 20,000rpm
Noise :	≤70dB
Spray :	External channel
Hole :	2 holes or 4 holes

TW-42 Teeth whitening Unit



Features:

1. Economic and cost-effective with simple design.
2. Dental chair mounted type and cabinet mounted type for choice.
3. LED blue light with minimized heat.
4. LCD display screen.
5. Wide range cover for whitening of 16 teeth.