

V2 Pro

Designed for High-End Clinics

Stable Structure

3-point ultra-stable chair frame structure.
Flexible backrest collision-resistant design.
5mm full-metal backrest + 12mm thickened base plate.
One-piece molded metal backrest.
Soft start and stop micro-sensation system.
Scientifically distributed force, low-noise operation.



Hygiene System

Microfiber antibacterial cushion, medical-grade stain-resistant.

Quick-release dual-face side box, magnetic touch design.

Imported medical-grade water and air tubing, resistant to high-frequency disinfection.

Removable ceramic cuspidor, high-quality full ceramic material.



Efficient Collaboration

673mm × 331mm Panoramic Wide-Range Workstation

** 12-key button control panel*

** 6-position hanging rotary arm system*

Multi-Function Side Control Tray

** 9-button design*

** Independent operation area*

Professional Lighting

8-piece imported LED dual-color dental light

** Adjustable color temperature from 4000K to 5500K*

** Even illumination of 8000–45000 LUX*

