

G1

The perfect solution for high demand

Firmness: All down to ground side box

Comfort: Grand & ergonomic designed seat

Excellence: Memory system

Profession: Twelve bulbs LED Sensor dual-colored lamp

Advance: Imported leather cushion



Floor standing side box, more stable and in line with implanting

synchronize seat frame(ergonomic comfort and zero back friction)

Stainless steel 304 tool tray, meeting the requirements of planting and daily disinfection

Trolley implant type

Optional Main Table



Luxury 20-LED dual color shadowless lamp



Top-mounted Type




Separate Cart Type

Seat Size	Large Width - 54cm , Seat - 64cm Backrest
Package	150*106*125cm , G.W : 290KG
Container Volume	10units/20GP , 32units/40HQ

Major Parts

1. Water heating system
2. Spittoon flushing and cup filling (supported by Italian solenoid valves)
3. Quiet electric lift chair (supported by Denmark motors)
4. LCD display control panel
5. Chair movement system with 9 memory positions
6. Panoramic LED film viewer
7. Assistant control panel
8. Rotatable ceramic spittoon
9. Shadowless LED sensor dental light
10. Multifunction pedal switch
11. Double pure water bottles, 600ml each
12. Glass tray and built-in tissue box
13. Deluxe dentist stool
14. Imported polyether inner tubing



Technical Specification

Supply voltage	: AC110V/220V 60Hz/50Hz
Motor voltage	: DC24V
Air pressure	: 0.5MPa - 0.8MPa
Water pressure	: 0.2MPa - 0.4MPa

Dentist Stool



01



02



03

Optional Dental Lamp



L1



L2



L3

Optional Accessory



Handpiece Set



Built-in Scaler



Built-in LED
Curing Light



Intra-oral
Camera Set



Medical
Compressor



Top-mounted
Main Table

Upholstery Color

Microfiber Upholstery



PU Upholstery

