

G5

Deluxe Experience


Firmness: All down to ground side box

Comfort: Three folding seat that flex at the same time offering more comfort

Humanization: Humanized functions for high effective work experience

Profession: 20 bulbs led surgical lamp

Safety: Multi-protection from the Interlock system



Luxury 20-LED dual color shadowless lamp

Optional Main Table



Top-mounted Type



Separate Cart Type

**HUMAN
CENTERED DESIGN**
PHILOSOPHY
CREATES HIGHER
VALUE



Three- folding seat system

Seat Size	Large Width - 54cm , Seat - 64cm Backrest
Package	150*106*125cm , G.W : 300KG
Container Volume	10units/20GP , 32units/40HQ

Major Parts
1. Water heating system
2. Spittoon flushing and cup filling (supported by Italian solenoid valves)
3. Quiet electric lift folding chair (supported by Denmark motors)
4. LCD screen control panel
5. Chair movement system with 9 memory positions
6. Panoramic LED film viewer
7. Assistant control panel
8. Rotatable ceramic spittoon
9. 8 LED surgical dental light
10. Multifunction pedal switch
11. Double pure water bottles, 600ml each
12. Glass tray and built-in tissue box
13. Deluxe dentist stool
14. Imported polyether inner tubing



Technical Specification
Supply voltage : AC110V/220V 60Hz/50Hz
Motor voltage : DC24V
Air pressure : 0.5MPa - 0.8MPa
Water pressure : 0.2MPa - 0.4MPa

Dentist Stool



01

02

03

L1

L2

L3

Optional Accessory



Handpiece Set

Built-in Scaler

Built-in LED Curing Light

Intra-oral Camera Set

Medical Compressor

Top-mounted Main Table

Upholstery Color

Microfiber Upholstery



PU Upholstery

