

Configuration ★ Standard / △ Optional

		Luxury	Standard			Luxury	Standard
Patient Chair	Seamless PU cushion	△	★	Main Side Box	90 degrees rotatable side box	★	★
	Luxury nappa leather cushion	★	△		Rotatable Ceramic Cuspidor	★	★
	Synchronized Chair & Chair Interlock System	★	★		Imported water & air tube	★	★
	Spitting / Treatment / Emergency Chair Position	★	★		Independent water filter, dual filtration function	★	★
	Soft-start and soft-stop function & backrest anti-collision	★	★	Dentist Stool	Luxury Dentist stool (MB1)	★	★
Instrument Tray	Instrument Tray heating Water System	★	△	Dental Light	Upgraded Assistant stool (MB3)	★	△
	9 sets of chair position memory	★	△		10 LED bead dental light	★	△
	6-position Adjustable Instrument Holder	★	★	Options	Reflector lamp (2 kinds of light mode)	△	★
	Handpiece anti-retraction function	△	△		Selected imported leather cushion	△	△
	One button Water-air-electricity on/off switch	★	★		Upgrade to mobile cart	△	△
	Handpiece independent water adjustment	★	★		Built-out floor box	△	△
	Built-out X-ray film viewer	★	★		High and low speed handpiece	△	△
Assustant Control	3-way Syringe (with warm water)	★	★		Built-in scaler and Built-in curing light	△	△
	9 keys control panel	★	★		Intra-oral camera and its holder	△	△
	Adjustable Strong/Weak Suction Head	★	★		Air Compressor	△	△
Floor Box	Built-in Floor Box	★	★				
	Current filter and Voltage display	★	★				



INNOVATIVE IDEA, EXQUISITE CRAFTSMANSHIP

We Medical is located in Nansha District, a key development area of Guangzhou City, with vast development potential and a solid industrial foundation. The company adheres to strict product inspection and process standards, with advanced production and quality control systems. As an emerging enterprise in dental chair manufacturing, MIPONT Medical is committed to innovation and development, providing high-quality and refined solutions for dental clinics.




DESIGNED FOR MODERN HIGH-END DENTAL CLINICS



TO CREATE A MODERN HIGH-END CLINIC

Designed for modern high-end dental clinics

Adhering to the design philosophy of "Less is More", We emphasizes the combination of elegant design and medical needs, infusing products with innovative designs tailored for modern high-end dental clinics, with strong insight of user needs and market needs.



We has established professional lab with a variety of lab equipment and over 30 sets of lab equipment, running over 150 tests covering more than 280 kinds of materials, including environmental testing and reliability testing and etc.

OVER 30 Professional lab equipment
OVER 50 Professional Quality Control team
OVER 150 Testing

Production and Management

We has a strict and lean production management system and a refined 5S management workshop to ensure the efficient operation of production, providing high quality product and dedicated service to our partners.

Strict Quality Control

We conducts steps of quality control from the raw material with product reliability test, simulation test of the use of dental chairs, extensive research and testing of the existing products, and product upgrades, to ensure well-designed and well-produced high-end product.

The high-end modern dental chair, more style of the clinic





M100 Luxury

Maintaining the design advantage as a leading position.

M100 Luxury is a high-end dental chair with good combination of high quality materials, good stability, comfort and durability, perfect choice for high-end dental clinic.

Cloud like comfort.


“ AS IF
LYING ON A SOFT CLOUD
”



COMFORTABLE TREATMENT FOR PATIENTS

Ultra-stable comfort soft activate system, soft start and stop function, Synchronized chair design, effectively reduce the patient's fear, let the patient experience immersive treatment; 380mm industry's lowest chair position, making it convenient for elderly and child patients.

Aluminum alloy armrests leather and cotton provide delicate touches.
Modern Nappa Leather cushion for extremely comfortable lying.
Curved headrest, better fitting the patient's head.



heating water function, provides warm water for handpiece and 3-way syringe.
Effective reduce the soreness during treatment.

Ultimate experience for dentist.

LUXURY DENTIST STOOL

We in-depth study of ergonomic design, comfortable lumbar support.
The stool's height and backrest recline are fully adjustable for optimal comfort.

Four-way adjustment, adjustable seat lift, recline back freely.
Silent pulleys, more durable aluminum alloy base frame.



PRECISE OPERATION FOOT PEDAL

Every press delivers accurate feedback

Low-position design for easy operation

Long durability

OUTSTANDING LIGHTING



LED HIGH-END SENSOR DENTAL LIGHT

Clearer visibility with 10 Philips-imported LED beads

Seamless adjustment | Excellent illumination | Precise control

Brightness memory function | 3 color modes | 7-80k Lux

Color temperature: 4000K-5700K | High CRI (>90Ra)

DUAL CONTROL DESIGN

SAFE AND EFFICIENT



FOR ASSISTANT

9-buttons multifunctional side control, providing efficient operation.
Comfortable nurse stool, adopting high-density molded cotton.

Strong metal backrest


3mm thickened full metal base for stronger support.

Movable side box & cuspidor

90 degrees rotate for easy operation.

Thicker curved plate, firmer connection

Extra large die-cast 30mm aluminum plate.



Center connection to the sidebox

Easier to install the equipment outward rotatable sidebox.
Provide flexibility for narrow space installation without splitting the sidebox.

Solid bottom plate

Thickness 12mm, stable as a rock. The unit weighs 250KG, making it more stable.



M100 Standard

Combine both quality and style.

KNOW YOU BETTER



FOR DENTIST

SELECTED REFLECTOR LAMP

Gentle light source reduces direct stimulation to patients' eyes. 100,000 times high-frequency life test. Imported Philips lamp LED.

Manual/sensing dual-mode for different treatment
2 modes: yellow and white light
8-30k Lux adjustable brightness



Large-sized instrument tray

Water-air-electricity with one on/off switch control.
Dust-proof and anti-fouling silicone pad.
Independent water adjustment.



Luxury dentist stool


Silent pulley, smooth movement.



FOR PATIENTS

SYNCHRONIZED CHAIR MOVEMENT FUNCTION

FOR MORE COMFORT DURING THE CHAIR ADJUSTMENT



LOWEST SEAT POSITION

380MM CHAIR POSITION
IDEAL FOR ELDERLY AND PEDIATRIC PATIENTS


CURVED DESIGN HEADREST

MORE SUITABLE FOR THE CURVATURE OF THE HUMAN HEAD STRUCTURE
OFFERING PATIENTS MORE COMFORTABLE AND SAFETY TREATMENT

FOR ASSISTANT

ASSISTANT INSTRUMENT TRAY

9-buttons multifunctional side control,
providing efficient operation.



Solid bottom plate

Thickness 12mm, stable as a rock. The unit weighs 250KG, making it more stable.

Center connection to the sidebox

Easier to install the equipment outward rotatable sidebox.
Provide flexibility for narrow space installation without splitting the sidebox.

Soft Start and Soft Stop

Chair movement without sudden quiver. More comfortable for patient.

Strong metal backrest

3MM thickened full metal base for stronger support.

Thicker curved plate

Extra large die-cast 30mm aluminum plate.





Long-term Solution For Clinic

Safety with anti-collision system

Emergency stop makes the operation safer when encountering obstacles.

Double water filtration

Double water filtration improves water quality and ensures water safety.

Imported water and air tube

Good flexibility, corrosion resisting, durable and high quality.

Durable Silent motor

Leading quality, quiet and stable.

Digital voltmeter

Display of input voltage number, additional filter, anti-voltage fluctuation.

Optimized design layout

Optimize the layout design of internal pipelines to effectively extend the service life of internal components.



M100 Luxury

UPGRADED

MOBILE CART

Elegant mobile cart with motor, Height adjustable.
Stainless steel instrument tray, durable and sterilizable.



Different Color
Different Style

Seamless PU Cushion (Standard for M100 Standard)



Nappa Leather (Standard for M100 Luxury)



Specially Selected Microfiber Leather (Optional)

