

S670 DENTAL UNIT

Simple and comfortable
Reliable to use



► LED oral light(Six LED beads with white light)

Induction on/off brightness adjustment to reduce eye fatigue and achieve accurate color comparison

Yellow light mode

Can manually switch yellow light to prevent resin from curing too quickly, providing more ample operating time;

White light mode

Adjustable color temperature, induction on/off brightness adjustment, meet the operational habits of different doctors;

Keep Improving Casting Quality

We are well aware of the importance of safety for dental chairs. Therefore, in the product design and manufacturing process, we strictly follow national and industry standards, use high-quality materials and advanced technology to ensure the safety and reliability of the products.



► Assistant rack combination device(Delayed shutdown)

Aluminum alloy handle detachable for cleaning and disinfection

► Assistant rack silicone pad

Anti scratch, easy to clean, and disinfectable

► Aluminum alloy instrument tray (interlocking)

The handle can be disassembled for disinfection, and the hanging rack can be rotated


► Instrument tray silicone pad

Anti scratch, easy to clean and disinfectable



► External water storage bottle device

One click conversion/switch, convenient and fast



► Detachable strong weak suction filtration device

Easy to clean



► Floor mounted box

Installed on the ground with high stability

S670

DENTAL UNIT

Engineering electric chair Comfortable seat cushion design

The seat design is based on the human body shape, providing perfect comfort and support, especially for patients undergoing long-term treatment.



► Silent electric chair (Protect)

Designed with a minimum seating position of 380mm, convenient for children and elderly



► Ceramic spittoon, automatic water supply system

Detachable/rotating, flushing and mouthwash quantitative water supply



► Emergency stop

One click stop, safe, and convenient



► Multifunctional foot pedal

Control up and down/supine/ cupping/ water supply/ oral light / spitting position/ mobile phone working water status/ speed



► Water/gas/electricity main switch device

One click on/off, convenient and fast (anti suction device)

Factory Configuration List		S610	S620	S630	S640	S640 Children chair	S670	S680 Planting chair
Electric chair	24V DC silent electric chair (safety protection)	■	■	■	■	■	■	■
	3 memory bit control systems (integrated circuit boards)	■						
	9 memory bit control systems (integrated circuit boards)		■	■	■	■	■	■
	Folding double joint headrest	■	■	■	■	■	■	■
	Dual armrest (rotatable/aluminum alloy material/engineering plastic cover)	■	■	■	■	■	■	■
	Water/gas/electricity main switch device (anti suction device)	■	■	■	■	■	■	■
	Seamed seat cushion (environmentally friendly leather)				■	■	■	■
	Seamless seat cushion (anti static PU leather)	■	■	■				
Instrument tray	Built-in floor box	■	■	■	■	■	■	■
	Hanging instrument tray (aluminum alloy seat/interlock/air lock device)			■				
	Underhung instrument tray (aluminum alloy seat/interlock/air lock device)	■	■				■	
	Cart type instrument tray (aluminum alloy seat/interlock/air lock device)				■	■		■
	Stainless steel tray (rotatable)			■	■	■		■
	Three purpose gun (cold/hot)	■	■	■	■	■	■	■
	Mobile phone hanging device (rotatable)	■	■	■	■	■	■	■
	LED film viewing light	■	■	■	■	■	■	■
Mouth lamp	Lamp arm device (with damping device)	■	■	■	■	■	■	■
	S1 LED oral light (anti resin curing mode/induction on/off/illumination adjustment)	■						
	S2 LED oral light (white light, yellow light, mixed light, induction on/off/illumination adjustment)		■	■		■		
	S3 LED oral light (white light, yellow light, mixed light, induction on/off/illumination adjustment)						■	
	S4 LED surgical light (white light, yellow light, mixed light, induction on/off/illumination adjustment)							■
Spittoon	Ceramic spittoon (rotatable/detachable)	■		■	■	■	■	■
	Glass spittoon cylinder (rotatable/detachable)		■					
	Automatic water supply system for quantitative flushing and mouthwash	■	■	■	■	■	■	■
Chassis	Assistant rack combination device (double joint foldable rotation)	■	■	■	■	■	■	■
	Multifunctional foot pedal (up and down/supine/sputum flushing/water supply/ light/mouth- wash/cooling water/blowing debris/ speed)		■	■	■	■	■	■
	Multifunctional foot pedal (up and down/supine/sputum flushing/water supply/ light/mouth- wash/cooling water/blowing debris)	■						
	Strong/weak suction device (detachable aluminum alloy handle/delayed closing)	■	■	■	■	■	■	■
	External water storage device	■	■	■	■	■	■	■
	Medical grade water/trachea	■	■	■	■	■	■	■
	USB connection cable (integrated circuit)	■	■	■	■	■	■	■
Doctor's chair	Constant temperature heating device	■	■	■	■	■	■	■
	S30 Doctor's Seat	■	■	■	■	■	■	■
	S60 Doctor's Seat		■	■	■	■	■	■
Choose	Waterway disinfection instrument	□	□	□	□	□	□	□
	External floor box	□	□	□	□	□	□	□
	Induction water supply device	□	□	□	□	□	□	□
	Multimedia oral endoscope (endoscope bracket)	□	□	□	□	□	□	□
	Central negative pressure device	□	□	□	□	□	□	□
	LED light estimation (built-in)	□	□	□	□	□	□	□
	Tooth cleaning machine (built-in)	□	□	□	□	□	□	□
	Imported/domestic high-speed mobile phones	□	□	□	□	□	□	□
	Imported/domestic low-speed mobile phones	□	□	□	□	□	□	□

Choose seat cushion colors

Seamed eco-friendly leather

H1

H2

H3

H4

H5

H6

H7

H8

H9

H10

Optional seat cushion colors

Seamless PU leather

P01

P03

P05

P06

P07

P09

P10

P11

Option



S1 LED oral light



S2 LED oral light



S3 LED oral light



S4 LED oral light



S4 LED surgical light (recording system)



Intraoral Camera system



Hanging instrument tray



Cart type instrument tray



S30 Doctor's Seat



S60 Doctor's Seat







UVC water sterilizer
Deep ultraviolet radiation
Overflow physical sterilization
Isolation without addition
More secure
Rapid sterilization rate:
Millisecond sterilization rate
99.99%



UVC Air disinfection
Deep ultraviolet radiation
Overflow physical sterilization
Isolation without addition
More secure
Rapid sterilization rate:
Millisecond sterilization rate
99.99% (Air)

Construction reference drawings





Technical Parameter

■ supply voltage	AC220V/50HZ
■ Water pressure of water source	0.2MPa-0.4MPa
■ Air Pressure	0.5MPa-0.8MPa
■ Air source displacement	>55L/min

page/30