






LED dental lamp

Yellow and white light,
Stepless induction regulation
And manual adjustment, adapt
Different treatment needs.



Multi-function foot control

Can be controlled by
chair lifting, backrest,
handpiece status and
spittoon flushing water
supply, improve the
doctor's work efficiency.

Multi-function foot control

Can be controlled by
chair lifting, backrest,
handpiece status and
spittoon flushing water
supply, improve the
doctor's work efficiency.



Top-mounted instrument tray

The treatment table is set
under the cushion, can
rotate left and right
in 170° range, more
convenient for doctors
to operate.





Aluminium alloy base

Precision casting process is adopted, made of aluminum alloy base, Steady, rust-proof, Anti-deformation, firm and stable.



Aluminium alloy back

The backrest is precision cast Made of aluminum alloy, the seat is more stable and firm. Suitable for all kinds of bodies.

Left hand operating system

Quick switching left and right hand operation System (optional), for different habits Doctors bring convenience.



Square assistant table

Wide and simple, Can rotate 120 °





**All tubes
from the USA**
Durable
Life for a long time
Environmental protection




**Handpiece tube
check valve**
Prevent sewage back
suction, avoid cross
infection of patients.

**Stainless steel
handpiece joint**
Silica gel handpiece tube,
acid and alkali resistant
Anti low temperature,
anti static electricity



**The abrupt
stop switch**
Stop in case of
emergency, prevent
medical errors





Independent air controlled water valve

- Each handpiece tube uses a separate air-controlled water valve, independent of each other.



Inlet pipe check valve

- The intake pipe adopts anti-suction device, Prevent sewage back suction, avoid pollution Urban water.

Pressure relief valve

- When the pressure is too large, it will be automatically relieved, which can protect the water bottle and make the work more safe.



Mute motor

- The motor runs smoothly and the noise is low. When the treatment is running, the feeling of pushing back is not strong. Reduce the patient's fear.

