

Packing list



Dentist Cart

- | | | |
|----|--------------------------------------|------------|
| 1 | Main Unit | *1 |
| 2 | User Manual | *1 |
| 3 | Power Cord | *1 |
| 4 | Water Tank | *1 |
| 5 | Guide Sleeve Plug | *1 |
| 6 | Motor Handpiece 1 (with cable) | *1 |
| 7 | Motor Handpiece 2 (with cable) | *1 |
| 8 | Three-way syringe (with cable) | *1 |
| 9 | Foot Pedal (with foot pedal hook) | *1 |
| 10 | Support Rod | *1 |
| 11 | Bracket | *1 |
| 12 | Power Fuse 6.3AH 250V | *2 |
| 13 | Certificate of Conformity | *1 |
| 14 | Detachable extension table | *2 |
| 15 | Ultrasonic Scaler Handpiece and Tips | (optional) |
| 16 | Air prophy (with cable) | (optional) |



Assistant Cart

- | | | |
|---|-------------------------------|----|
| 1 | Main Unit | *1 |
| 2 | User Manual | *1 |
| 3 | Power Cord | *1 |
| 4 | High-performance suction hose | *2 |
| 5 | Saliva suction hose | *1 |
| 6 | Power Fuse 6.3AH 250V | *2 |
| 7 | Certificate of Conformity | *1 |

MOBILE DENTAL SOLUTION
NO NEED AIR COMPRESSOR AND PLUMBING

Implant & surgery | General practice (1:5 & 1:1) | Air polisher | Scaler | Three — way syringe



Dentist cart


Multi functions in one device

(Implant, restoration, periodontal, surgery, ect.)

- Electric driven
- Screen touch operation with digital display
- Convenient for mobile treatment
- Large tray design for instrument placement
- Silent (≤60dB)



Assistant cart



Clinical pictures



Colour available for dental chair



Clinical pictures

(End-user experience)



Application Scenario




Laying pipelines Wasting Time, Labor, Money



Plug and play (only electricity)



Much space occupation by different devices



Multiple uses in one machine




Fixed assembly in clinic



Portable/Mobile Treatment



Limited control area



Flexible left or right handed operation

Model	DS70
Size (excluding water tank)	650mm (L) × 550mm (W) × 910mm (H)
Net weight	45 KG
Input voltage	110/220VAC
Rated current	0.1A-4A
Power fuse	2*250V-T6.3AH
Maximum power consumption	400VA
Water tank size	400mm (L) × 140mm (W) × 220mm (H)
Fresh water tank capacity	5L
Cooling water flow at 100%	50 mL/min or more
Recommended water tank coolant	distilled water/purified water
Maximum allowable suspension load	≤1.5KG



Model	Suction 10
Size (excluding water tank)	720mm(L)×270mm(W)×820mm(H)
Net weight	48.5Kg
Input voltage	110/220VAC
Rated current	0.1A-5A
Power fuse	2*250V-T6.3AH
Maximum power consumption	500VA
Sewage tank size	350mm(L)×170mm(W)×280mm(H)
Sewage tank capacity	10L(Automatic overflow protection)
Suction pump	Piston Pump
Ultimate Negative Pressure Value	≥0.06Mpa
Pumping Rate	Strong Suction: 10L/min, Weak Suction: 4.5L/min