

15 years
estimated service lifetime

Long-lasting quality
Extraordinary comes from every detail



Coexistence of looks

Grasp the Simplicity of Modern Clinic

Streamlined Design
Depicting Modern Beauty

V2 Configuration

Patient Chair	Comfortable skin-friendly leather																									
	Fashionable double armrests																									
	Spitting reset, treatment chair position																									
	High-quality domestic motor																									
	Adjustable headrest																									
	Multifunctional foot pedal																									
Instrument Tray	6-position rotatable handpiece holder																									
	14-key functional control panel																									
	One-button water, air and electricity switch																									
	External X-ray film viewer																									
	Round pressure gauge																									
	Original blister cover																									
Side Box	Stainless steel handle																									
	Double-door side box																									
	Rotating and removable ceramic spittoon																									
Assistant Tray	High-quality water air tube	Optional seat cushion colors <table><tr><td>V01</td><td>V02</td><td>V03</td><td>V04</td></tr><tr><td></td><td></td><td></td><td></td></tr><tr><td>Sky blue</td><td>Light blue</td><td>Medium blue</td><td>Light gray</td></tr><tr><td>V05</td><td>V06</td><td>V07</td><td>V08</td></tr><tr><td></td><td></td><td></td><td></td></tr><tr><td>Greenish brown</td><td>Light red</td><td>Light orange</td><td>Apple green</td></tr></table>	V01	V02	V03	V04					Sky blue	Light blue	Medium blue	Light gray	V05	V06	V07	V08					Greenish brown	Light red	Light orange	Apple green
	V01		V02	V03	V04																					
																										
Sky blue	Light blue		Medium blue	Light gray																						
V05	V06		V07	V08																						
																										
Greenish brown	Light red		Light orange	Apple green																						
3-way syringe (warm)																										
9-key control																										
Adjustable strong/weak suction																										
Operation Light	6-bulb LED light (yellow and white color with induction)																									
Floor Box	External floor box																									
Dentist Stool	Advanced multi-directional adjustable dentist stool (skin-friendly leather)																									

Creating High-End Clinic Experience

Professional dental equipment manufacturer, innovating modern dental chair



Pursuit for Perfection

We are keen on every details to be perfect

Rotating spittoon

Closer and easier for patients
Deeper spittoon with a unique design, effectively prevent-
ing water splashing
Detachable glass spittoon, easier for disinfection
(Ceramic is optional)

Luxury dentist stool

Adjustable in four directions

Full of security

Chair interlocking function, more safe during treatment

Rotatable bracket

user-friendly design
Flexible and Adjustable according to personal habits

Easy to clean

Visible External Suction Filter
Convenient for daily cleaning

External ground box

Flexible installation,Easy to clean



Comfort Enhanced

Exquisitely comfortable for every step

Skin-friendly Leather Cushion

Comfortable, soft, and wear-resistant
Delicate touch, wear-resistant and scratch-resistant
Soft and breathable surface
comfortable for long-time treatment

Exquisite and comfortable

Curved headrest with adjustable
angles for easy positioning

Comfortable backrest

Pear-shape backrest perfectly fits patient's back

PU armrest

Exquisite and smooth

Practicality First

Every operation is effortless

Multi-functional Instrument Tray

Large-sized tray, providing more storage space

13 Keys Multifunction Control Panel

Control for dental lamp, treatment position
spit position, and Reset

Bottom mounted mechanical disc

Hanging tray providing more efficiency

- One-Key Water, Air, and Power Switch
- 6-Position Rotatable Brackets
- External X-ray film viewer

Exclusive Design Dental Lamp

Sensor and touch control,Superior lighting effect for various dental operation


- 6-Bulb Imported LED Light Source
- 3 Types of Light Source Modes
- Ra>90



Self-developed multi-functional foot pedal

A composite foot pedal with sputum collection
position, oral lamp, suction, and water supply function keys.

Comfortable control, no fatigue even after long-time use,
optimized pedal control angle and feedback force to
improve doctor efficiency.



Multi functional assistant control

8-key function control
Four handed operation is more efficient

Pneumatic control
Strong and weak suction handle with silicone

